


## Applications

Selective hydrolysis of amides in aqueous solution



## Kit description

The kit contains 6 diverse pre-formulated Amidase (AMD) biocatalysts as lyophilised powders, as well as pre-prepared phosphate buffer.

### AMDs contained in the screening kit

AMD 101
AMD 106
AMD 107
AMD 109
AMD 110
AMD 125

### Contents

AMDs	6 enzymes (50 mg)
Dithiothreitol (DTT)	100 mg
0.1M KH <sub>2</sub> PO <sub>4</sub> buffer (pH 7.2)	60 mL

## Screening Procedure


1. Label 6 x 2mL vials corresponding to the six different AMDs provided in the kit.
2. Add DTT (15.4 mg) in 10 mL of buffer provided to make a 1 mM DTT buffer\*\*.
3. Add 10 mg of enzyme into the labelled vial.
4. Add 450 µL of buffer containing 1 mM DTT to each vial.
5. Add 2-10 mg of amide in 50 µL of DMSO to each vial.
6. Incubate overnight at 37 °C with agitation.
7. The reaction can be analysed by HPLC or TLC to check for hydrolysis.

\*\*It is recommended to make the reaction mix solution fresh and use immediately. Avoid storage of the reaction mix as a solution, as this will degrade over time. An adequate supply of DTT and buffer is provided for five screens. Additional DTT or buffer can be purchased from Almac if required.

**Storage:** Recommend refrigeration at 4°C to preserve enzyme activity

### selectAZyme Offerings

- An ever-expanding biocatalysis team including molecular and microbiologists, enzymologists, bioinformaticians, organic chemists and analysts, all equipped with state-of-the art facilities.
- Expertise in gene identification, expression, fermentation and enzyme production, followed by the efficient use of enzymes to produce complex chiral APIs.
- Enzyme evolution based on computational re-design, semi-rational and random mutagenesis approaches, allowing access to bespoke biocatalysts with enhanced activity, selectivity and process robustness.
- Fully integrated biocatalyst development through screening, (chemo-) enzymatic route definition, process development and scale up (pilot plant facilities available).
- Rapid implementation of enzymatic steps in complex, multi-stage syntheses, leading to significant improvements in production yields and timelines.
- A simple business model that avoids IP issues.



### The selectAZyme Range of Enzyme Screening Kits

Our selectAZyme kits include a detailed user guide and come with all buffers, cofactors, recycling systems and reagents necessary to perform screens using standard laboratory equipment.

#### Carbonyl Reductase (CRED) biocatalysts

96 CRED biocatalysts for the production of chiral alcohols and/or use in cofactor recycling schemes

#### Aldehyde Reductase (ARED) biocatalysts

16 ARED biocatalysts

#### Hydrolase biocatalysts

48 commercially available hydrolases for selective acylation of alcohols and amines.

#### Nitrilase and Nitrile Hydratase (NHase) biocatalysts

9 NHases and 15 nitrilases

#### Transaminase (TAm) biocatalysts

96 TAm for the production of chiral amines from pro-chiral ketones.

#### Ene Reductase (ERED) biocatalysts

143 ERED biocatalysts for asymmetric reduction of activated alkenes

#### P450 Monooxygenase biocatalysts

96 P450 monooxygenase biocatalysts for a huge range of highly selective oxidations

### Want Almac to do the screening for you?

- Our experienced biocatalysis team can screen all of our enzymes against your target substrate(s) and simply provide the results.
- Flexible options for subsequent enzyme supply, evolution services, process development and scale up as required.

### Technical Contacts:

Prof. Tom Moody, Tel: +44 (0)28 3833 2200 Ext. 5517, E-mail: [tom.moody@almacgroup.com](mailto:tom.moody@almacgroup.com).

Dr. Derek Quinn, Tel: +44 (0)28 3833 2200 Ext. 5833, E-mail: [derek.quinn@almacgroup.com](mailto:derek.quinn@almacgroup.com).

Address: Almac Biocatalysis & Isotope Chemistry Group,

20 Seagoe Industrial Estate, Craigavon BT63 5QD

Web: [www.almacgroup.com](http://www.almacgroup.com),

Email: [biocatalysis@almacgroup.com](mailto:biocatalysis@almacgroup.com)