


### Applications

Deoxyribose-5-phosphate Aldolases (DERAs) catalyse a reversible aldol reaction between acetaldehyde and D-glyceraldehyde 3-phosphate to generate 2-deoxy-D-ribose 5-phosphate.



### Kit description

The kit contains 96 diverse pre-formulated DERAs biocatalysts as lyophilised powders, as well as pre-prepared Tris buffer.

#### DERA contained in the screening kit

	1	2	3	4	5	6	7	8	9	10	11	12
A	1	9	17	25	33	41	49	57	65	73	81	89
B	2	10	18	26	34	42	50	58	66	74	82	90
C	3	11	19	27	35	43	51	59	67	75	83	91
D	4	12	20	28	36	44	52	60	68	76	84	92
E	5	13	21	29	37	45	53	61	69	77	85	93
F	6	14	22	30	38	46	54	62	70	78	86	94
G	7	15	23	31	39	47	55	63	71	79	87	95
H	8	16	24	32	40	48	56	64	72	80	88	96

#### Contents

DERAs	96 enzymes (10 mg each)
DMSO	1 vial (10 mL)
0.1M Tris buffer (pH7.4)	1 bottle (50 mL)

### Screening Procedure

1. Add 500  $\mu$ L Tris-HCl buffer to each well.
2. Add a solution of  $\sim$ 5 mg substrate and 2 equivalent of acetaldehyde or acetone. If needed it can be dissolved in organic solvent (50-150  $\mu$ L, depending on solubility) such as DMSO or MTBE.
3. Shake/stir at room temperature (or ideally 35  $^{\circ}$ C). Agitate overnight.
4. Extract product with an organic solvent (MTBE, EtOAc etc.), if needed.
5. Analyse sample by GC/HPLC to determine conversion.

\*It is recommended to make the reaction mix solution fresh and use immediately. Avoid storage of the reaction mix as a solution, as this will degrade over time.

**Storage:** Recommend refrigeration at 4 $^{\circ}$ C to preserve enzyme activity.




# Aldolase (DERA)

## Enzyme Screening Kit

### DERAESK-9600

### selectAZyme Offerings

- An ever-expanding biocatalysis team including molecular and microbiologists, enzymologists, bioinformaticians, organic chemists and analysts, all equipped with state-of-the art facilities.
- Expertise in gene identification, expression, fermentation and enzyme production, followed by the efficient use of enzymes to produce complex chiral APIs.
- Enzyme evolution based on computational re-design, semi-rational and random mutagenesis approaches, allowing access to bespoke biocatalysts with enhanced activity, selectivity and process robustness.
- Fully integrated biocatalyst development through screening, (chemo-) enzymatic route definition, process development and scale up (pilot plant facilities available).
- Rapid implementation of enzymatic steps in complex, multi-stage syntheses, leading to significant improvements in production yields and timelines.
- A simple business model that avoids IP issues.



### The selectAZyme Range of Enzyme Screening Kits

Our selectAZyme kits include a detailed user guide and come with all buffers, cofactors, recycling systems and reagents necessary to perform screens using standard laboratory equipment.

#### Carbonyl Reductase (CRED) biocatalysts

96 CRED biocatalysts for the production of chiral alcohols and/or use in cofactor recycling schemes

#### Aldehyde Reductase (ARED) biocatalysts

16 ARED biocatalysts

#### Hydrolase biocatalysts

48 commercially available hydrolases for selective acylation of alcohols and amines.

#### Nitrilase and Nitrile Hydratase (NHase) biocatalysts

9 NHases and 15 nitrilases

#### Transaminase (TAm) biocatalysts

96 TAm for the production of chiral amines from pro-chiral ketones.

#### Ene Reductase (ERED) biocatalysts

143 ERED biocatalysts for asymmetric reduction of activated alkenes

#### P450 Monooxygenase biocatalysts

96 P450 monooxygenase biocatalysts for a huge range of highly selective oxidations

### Want Almac to do the screening for you?

- Our experienced biocatalysis team can screen all of our enzymes against your target substrate(s) and simply provide the results.
- Flexible options for subsequent enzyme supply, evolution services, process development and scale up as required.

### Technical Contacts:

Prof. Tom Moody, Tel: +44 (0)28 3833 2200 Ext. 5517, E-mail: [tom.moody@almacgroup.com](mailto:tom.moody@almacgroup.com).

Dr. Derek Quinn, Tel: +44 (0)28 3833 2200 Ext. 5833, E-mail: [derek.quinn@almacgroup.com](mailto:derek.quinn@almacgroup.com).

Address: Almac Biocatalysis & Isotope Chemistry Group,

20 Seagoe Industrial Estate, Craigavon BT63 5QD

Web: [www.almacgroup.com](http://www.almacgroup.com),

Email: [biocatalysis@almacgroup.com](mailto:biocatalysis@almacgroup.com)