



## Applications

Selective hydrolysis of phosphotriester compounds in aqueous media.



## Substrate range

A wide variety of structurally diverse phosphotriesters. A selection is shown below.



## Kit description

The kit contains 8 diverse pre-formulated phosphotriesterase (PTE) biocatalysts as lyophilised powders.

### PTEs present in enzyme kit

PTE-001
PTE-002
PTE-003
PTE-004
PTE-005
PTE-006
PTE-007
PTE-008

### Contents

PTEs	8 enzymes (50 mg each)
0.1M Tris buffer (pH 7.4)	1 bottle (50 mL)

## Screening Procedure

1. Add into a flask/vial, add 10 mg phosphotriesterase enzyme in 500  $\mu\text{L}$  Tris-HCl buffer pH 7.4\*.
2. Add a solution of  $\sim 1$  mg substrate. If needed, the substrate can be dissolved in organic solvent such as DMSO or MTBE (5-10% of the final volume).
3. Add buffer to a final concentration of 1 mL.
4. Shake/stir at room temperature (or ideally 30  $^{\circ}\text{C}$ ). Agitate overnight.
5. Extract product with an organic solvent (MTBE, EtOAc etc.), if needed.
6. Analyse sample by GC/HPLC to determine conversion.


\*It is recommended to make the components fresh and use immediately.

**Storage:** Recommend refrigeration at 4 $^{\circ}\text{C}$  to preserve enzyme activity.



## selectAZyme Offerings

- An ever-expanding biocatalysis team including molecular and microbiologists, enzymologists, bioinformaticians, organic chemists and analysts, all equipped with state-of-the art facilities.
- Expertise in gene identification, expression, fermentation and enzyme production, followed by the efficient use of enzymes to produce complex chiral APIs.
- Enzyme evolution based on computational re-design, semi-rational and random mutagenesis approaches, allowing access to bespoke biocatalysts with enhanced activity, selectivity and process robustness.
- Fully integrated biocatalyst development through screening, (chemo-) enzymatic route definition, process development and scale up (pilot plant facilities available).
- Rapid implementation of enzymatic steps in complex, multi-stage syntheses, leading to significant improvements in production yields and timelines.
- A simple business model that avoids IP issues.



## The selectAZyme Range of Enzyme Screening Kits

Our selectAZyme kits include a detailed user guide and come with all buffers, cofactors, recycling systems and reagents necessary to perform screens using standard laboratory equipment.

### Carbonyl Reductase (CRED) biocatalysts

96 CRED biocatalysts for the production of chiral alcohols and/or use in cofactor recycling schemes

### Aldehyde Reductase (ARED) biocatalysts

16 ARED biocatalysts

### Hydrolase biocatalysts

48 commercially available hydrolases for selective acylation of alcohols and amines.

### Nitrilase and Nitrile Hydratase (NHase) biocatalysts

9 NHases and 15 nitrilases

### Transaminase (TAm) biocatalysts

96 TAm for the production of chiral amines from pro-chiral ketones.

### Ene Reductase (ERED) biocatalysts

143 ERED biocatalysts for asymmetric reduction of activated alkenes

### P450 Monooxygenase biocatalysts

96 P450 monooxygenase biocatalysts for a huge range of highly selective oxidations

## Want Almac to do the screening for you?

- Our experienced biocatalysis team can screen all of our enzymes against your target substrate(s) and simply provide the results.
- Flexible options for subsequent enzyme supply, evolution services, process development and scale up as required.

## Technical Contacts:

Prof. Tom Moody, Tel: +44 (0)28 3833 2200 Ext. 5517, E-mail: [tom.moody@almacgroup.com](mailto:tom.moody@almacgroup.com).

Dr. Derek Quinn, Tel: +44 (0)28 3833 2200 Ext. 5833, E-mail: [derek.quinn@almacgroup.com](mailto:derek.quinn@almacgroup.com).

Address: Almac Biocatalysis & Isotope Chemistry Group,

20 Seagoe Industrial Estate, Craigavon BT63 5QD

Web: [www.almacgroup.com](http://www.almacgroup.com),

Email: [biocatalysis@almacgroup.com](mailto:biocatalysis@almacgroup.com)