

CASE STUDY

Trace Impurity Identification using Q-TOF Mass Spectroscopy

Authors: Dr. Osama Chahrour, Alan Thompson and Dr. John Malone

Accurate Mass Spectrometry (MS) techniques are ideally suited to aid in the structural elucidation of trace impurities. A rapid approach used for the isolation and characterisation of trace impurities in drug substances and formulated drug products has been developed. The approach has been successfully employed for the characterisation of impurities from different substances covering a wide variety of structural classes.




BACKGROUND

FDA and International Conference on Harmonization (ICH) set strict requirements regarding impurities in APIs and formulated drug products. It defines thresholds at which impurities must be identified or adequately tested. Investigations are initiated during the early stages of drug discovery or product development such as the pre-investigational new drug (pre-IND) stage.

Accurate Mass Spectrometry (MS) techniques are ideally suited to aid in the structural elucidation of trace impurities. A rapid approach used for the isolation and characterisation of trace impurities in drug substances and formulated drug products has been developed. The approach has been successfully employed for the characterisation of impurities from different substances covering a wide variety of structural classes. A case study is presented here in.

Outline of the impurity isolation and structure elucidation work flow


Molecular Feature Extraction (MFE) and molecular formula generation (MFG) algorithms, along with MassHunter Molecular Structure Correlator (MSC) greatly aid impurity identification and structure elucidation using accurate mass and targeted MS/MS data.



Impurity ID by accurate mass LC-MS/MS (LC-ESI-QTOF)


Almac completed a laboratory feasibility study and manufactured a demonstration batch of API. The work involved necessary analytical and process development to enable safe and reliable process scale up and API stability evaluation conforming to cGMP.

Structure of API




Optimised accurate mass impurity profiling successfully identified the molecular formula of all impurity peaks present in the API chromatography with high score using Molecular Feature Extraction (MFE) and Molecular Formula Generation (MFG) algorithms.

Accurate mass impurity profiling




Almac conducted a stability trial on the API to investigate a range of diluents, stabilisers and conditions. Impurity A concentration increased upon incubation and structural elucidation was necessary.

Stability study revealed increases in Impurity A




MS compatible method was developed and total ion chromatogram was collected (ESI-).

Total ion chromatogram



Blank subtraction and peak extraction of the previous chromatogram yielded the following accurate mass spectrum

Mass spectrum of impurity A



The molecular formula generation (MFG) algorithm successfully identified the formula with high score reflecting accurate-theoretical mass matching and the excellent isotopic distribution matching.


Compound Table

| Compound Label | RT | Mass | Formula | MFG Formula | MFG Diff (ppm) | DB Formula |
|-------------------|--------|---------|------------|-------------|----------------|------------|
| Cpd 1: C40 H66 O5 | 12.499 | 626.493 | C40 H66 O5 | C40 H66 O5 | -3.16 | C40 H66 O5 |

Predicted Isotope Match Table


| Isotope | m/z | Calc m/z | Diff (ppm) | Abund % | Calc Abund % | Abund Sum % | Calc Abund Sum % |
|---------|----------|----------|------------|---------|--------------|-------------|------------------|
| 1 | 625.4857 | 625.4837 | -3.09 | 100 | 100 | 63.86 | 63.86 |
| 2 | 626.4896 | 626.4872 | -3.96 | 44.65 | 44.2 | 28.51 | 28.23 |
| 3 | 627.4920 | 627.4903 | -2.73 | 10.28 | 10.56 | 6.57 | 6.75 |
| 4 | 628.4941 | 628.4933 | -1.15 | 1.65 | 1.79 | 1.06 | 1.14 |

Predicted isotope match



MS/MS fragmentation pattern @ 30 eV displayed three major daughter ions. The Molecular Formula Generation (MFG) algorithm successfully identified the daughter ions formula enabling the elucidation of the structure of the impurity.


Targeted MS/MS fragmentation



Peak List

| m/z | z | Abund | Formula |
|----------|---|---------|------------|
| 221.1552 | 1 | 51347.2 | C14 H21 O2 |
| 222.1586 | 1 | 8394 | |
| 303.2339 | 1 | 87005.1 | C20 H31 O2 |
| 304.2371 | 1 | 21094.2 | |
| 321.2444 | 1 | 65634.5 | C20 H33 O3 |
| 322.2475 | 1 | 15864.4 | |

Elucidation of Impurity A structure



Conclusion

The structure of the trace impurity A was successfully determined using Accurate Mass analysis by a Q-TOF mass spectrometer. The targeted MS/MS algorithm successfully identified the general chemical formulas of the impurity molecular and daughter ions. The information was used to elucidate the detailed chemical structure.

All our clients have unique needs.
That's why we develop unique solutions.

This is the ALMAC TOUCH™



GET IN TOUCH

UK

Almac Group
(Global Headquarters)
20 Seagoe Industrial Estate
Craigavon
BT63 5QD
United Kingdom

sciences@almacgroup.com
+44 28 3833 2200

US

Almac Group
(US Headquarters)
25 Fretz Road
Souderton, PA 18964
United States of America

sciences@almacgroup.com
+1 215 660 8500