



**CHEMTRIX**  
Scalable Flow Chemistry

## APPLICATION OF CONTINUOUS FLOW

Controlling Hazardous Processes from R&D to Production



Dr Charlotte Wiles – June 2023

## How are Flow Reactions Performed?

Solutions (typically) of reagents are pumped into a reactor, where they are;

- Mixed
- Heated or cooled
- Reacted for a specified period of time
- Collected for analysis or product isolation

Conceptually, a basic flow reactor comprises of;



Key is to understand the requirements of your process & design the reaction set-up accordingly

- This contrasts with batch practices where you typically adapt the process to suit the available vessel(s)!


## When to Consider Flow Chemistry?

A significant aspect of flow chemistry relates to performance versus scale. Reducing the reactor size can have a beneficial impact on;

- Heat transfer
- Mixing speed & mixing shear
- Repeatability
- Safety

Addressing priority areas of;

- Reduced energy costs & waste generation  
→ *Process sustainability*
- Increased reaction efficiency & process safety  
→ *Process intensification*



You must understand how these parameters influence your process before attempting scaling-up!

## Advanced Process Control (APC)

Compared to batch reactions, whereby changes evolve over hours, continuous processes can have reaction times in sec's to min's & operate under a state of control

- Analytical feedback therefore needs to be faster for continuous processes than batch
- Techniques that accurately reflect the process state need to be used
- After start-up, in flow you are typically monitoring for no process change

### 1. Design a manufacturing process to meet target(s)


- Identify & control critical quality attributes
- Understand impact of process parameter variability
- Select hardware to ensure target conditions are met

### 2. Monitor the process parameters to ensure consistent quality


- Dramatic reduction in operating costs
- Reduced reliance on man-power for plant operation

Replication of the unit & controls = same output quality


- Enabling decentralised manufacturing to be realised



## Process Control - Operating Under a State of Control



Many different PAT tools are used, ranging from T / P / flow / conductivity to FTIR, Raman, UV, MS & NMR




- If you only analyse the collected material you can miss these disturbances!

## Small, Agile Production Units

Flexibility is often requested when looking at small-scale manufacturing units. Flexibility can be towards;



- Raw material & product type, production location & production capacity



By adopting a modular approach, small flow production units can accommodate these needs & be replicated at different sites

- Whilst there is no 'one-size fits all' flow process
  - Similar reaction types can use a common platform (dosing, reactor etc.)
- Build in redundancy to enable a faster response to surges in material demand
  - Rapid, local ramp-up of manufacturing output

## From Lab to Manufacturing



### 1. Theoretical Evaluation

- Definition of process targets
- Identification of the 'problem'
- Theoretical opportunities for process improvement

### 2. Basic Feasibility Study

- Lab-scale trial(s)
- Deliver proof of concept
- Outline milestones for process optimisation


### 3. Parameter Optimisation

- Evaluate suitable hardware
- Demonstrate hardware feasibility
- Define process parameters for engineering design

### 4. Realisation

- Design pilot / production plant
- Build & test
- Commission

## A Safe & Efficient Method for the Preparation of Organic Azides




- No headspace – no  $\text{HN}_3$  build-up
- Safe access to 'Novel Operating Windows'
- Short reaction times, no decomposition

Demonstrating the rapid development of continuous flow protocols, using mg's of materials

Labtrix® Start – glass flow reactor(s)



- -20 to +195 °C at 20 bar







## A Safe & Efficient Method for the Preparation of Organic Azides

Advantages:

- Low consumption of material (mg's)
- Fast data generation
- Insight into key influences on product / by-product formations
- Early indication of viability of manufacture



## Customer Case – Biphasic, Dakin Oxidation




Drivers for flow method development;

- Safe use of  $\text{H}_2\text{O}_2$  at scale
- Shorter synthetic route
- Reduced reagent use
- Increased production rate

Protrix® Flow Reactor:

- Metal & glass-free reactor
- Integrated thermal control

→ 30 s reaction cf. 6 h in batch affording quant. catechol ( $4.8 \text{ kg } 8 \text{ h}^{-1}$ )




Partners:



## Customer Case - Wolff-Kishner Reduction

CDK9 Inhibitor, potential treatment of haematological malignancies (in Phase 1 clinical trials)



### Reasons for Flow:

- Safety risk as scale – hot N<sub>2</sub>H<sub>4</sub> in excess
- Challenging corrosivity
- Increased throughput / decreased reaction time

### Protrix<sup>®</sup> for Development:

- Metal & glass-free reactor with integrated thermal control
- Reduced N<sub>2</sub>H<sub>4</sub> excess due to no losses (3.5 eq. reduction)
- Time reduced from 2 h to 20 min
- Yield increased by 40 % (80 % with >99 % purity)



CCFLOW Partners: **GOFLOW**  **AstraZeneca** 



## Customer Case – Securing a Key Intermediate used in Vitamin B<sub>1</sub> Synthesis

3-Chloro-4-oxopentyl acetate is a key intermediate used in the synthetic preparation of Vitamin B<sub>1</sub> (Thiamine)

→ Vitamin B<sub>1</sub> is an antineuritic vitamin used in humans & animals

Drivers for flow process development;

- Improved yield of the chlorination step
- Reduced isolation costs
- Reduced waste generation



Using Protrix® Prof. Chen demonstrated the synthesis of 3-acetyl-3-chlorodihydrofuran-2(3H)-one;

- Reporting a 93 % yield (98 % purity), with a reaction time of 30 s (at 25 °C)
- Subsequent decarboxylative acylation afforded the target in 90 % yield (96 % purity)


Partner:



## Customer Case – Local Manufacture of Diabetes API's

Diabetes mellitus (DM) is a global health concern, with Type 2 DM most prevalent & managed with oral anti-diabetic drugs

→ Incidences of DM are growing in low- & middle-income countries resulting in increasing costs for its management



Drivers for flow method development;

- Control of API costs by manufacturing locally (currently API formulation only)
- Shorter, more efficient synthetic routes
- Increased production rate within a small footprint





Partner: **NELSON MANDELA**  
UNIVERSITY

## Customer Case – Local Manufacture of Diabetes API's

Using KiloFlow® glass flow reactors, Researchers at NMU (ZA) demonstrated the development & scale-up of four DM API's

→ Realising productivities ranging from 3.5 to 31.5 kg/week from a laboratory fume cupboard



Anti-diabetic Drug (API)	R Group	Day <sup>-1</sup>
 <p>R = (a) Glipizide; (b) Glibenclamide</p>	<p>(a) </p> <p>(b) </p>	0.6 kg
 <p>R = (c) Tolubutamide; (d) Chlorproamide</p>	<p>(c) CH<sub>3</sub></p> <p>(d) Cl</p>	4.4 kg
		4.5 kg

- Significant steps towards local API manufacture in South Africa!

Partner: **NELSON MANDELA**  
UNIVERSITY

## Customer Case – Manufacturing using a Lachrymator

Buchem BV identified an opportunity to improve their productivity for an existing product - combining batch & flow

Heart of the system was a 170 ml Plantrix® flow reactor containing SiC reactors, selected for chemical compatibility towards the challenging process & continuous mixing required, a 20 L Buchi rotavapor & a 50 L glass reactor followed

Cooperation: Buchem, Flowid & Chemtrix

Buchem – process & chemistry

Flowid – system design & engineering

Chemtrix – proof of concept & reactor design

Advantages:

- Higher productivity & robust process
- Smaller equipment
- Reduced material inventory
- Improved process control

→ Annual production requirements were met from a fume hood



Partners:  **BUChem BV** 

## Customer Case – Multi-step API & Intermediate Production

Fast process scale-up from R&D to production drove the innovation in this project!

Realised by combining Plantrix® flow reactors & an Agitated Tube Reactor (ATR) as a reactor/extractor. Smooth operation was ensured using super metering pumps.

Cooperation: Cipla, Pi Process Intensification, AM Technology, Fuji Techno & Chemtrix

### Result:

- Multi-product, multi-step flow chemistry production plant
- Compact, retrofitted into an existing asset base

### Advantages:

- Process flexibility & safety
- Speed of scale-up & enhanced process control via automation
- Reduced inventory of hazardous reagents
- Reduced 'batch to batch' variation



Partners:

**Cipla**



 **AMTechnology**

  
Fuji Techno Industries Corporation

 **processintensification**



## Customer Case – Selective Nitration



- Patheon uses flow reactors made of 3M™ SiC in a pharmaceutical production plant



Partner: **patheon**  
by Thermo Fisher Scientific

## Customer Case – Merck Sharpe Dohme (Gefapixant)



Standalone synthetic unit with monitoring of CPP utilising pressure, temperature & spectroscopy

- Operating at litres/min



© Copyright – Flowid

- Consistent product quality from continuous systems reduces the reliance on batch release & therefore work in progress

## Customer Case – Ajinomoto Bio-Pharma Services (CDMO)

Targeting a reduction in unit operations & increased productivity, Ajinomoto Bio-Pharma Services (BE) switched to flow!

### Customer Process:

- Three exothermic steps, requiring slow reagent addition (8 h)
- Potent product (OEB-5; 0.1-1  $\mu\text{g}/\text{m}^3$ )


Developing a multi-stage continuous process, at the end of 2017, gave:

- Increased process safety due to thermal control
- Minimised operator exposure
- Higher productivity *cf.* batch (x2)

### Technical Advantages:

- Small, mobile equipment tested in the lab & moved to production
- Reduced material inventory
- Enhanced process robustness *via* automation

→ Converting 70 % of the process to continuous improved batch capacity utilisation



Partner:



**BIO-PHARMA**  
SERVICES

## Customer Case – Ajinomoto Bio-Pharma Services (CDMO)

Targeting a reduction in unit operations & increased productivity, Ajinomoto Bio-Pharma Services (BE) switched to flow!



Scaled-up process commissioned 2020:

- Three chemical steps (2x continuous & 1x batch)
- >500 MT of active compound produced in the first year!

In addition to the advantages realised at the pilot-scale, the Industrial multi-stage continuous process enabled:

- Target price to be reached
- Increased annual output *cf.* dedicated synthesis unit (250 MT/y)
- Release of two 6000 L vessels for additional batch processes

Outlook:

- Downstream batch improvements increase output to >500 MT/y
- Development work shows 3<sup>rd</sup> step is also feasible in flow

Next Steps: Additional installation on a 2<sup>nd</sup> ABPS site in 2023 with all three steps continuous!

Partner:



**BIO·PHARMA**  
SERVICES

## Customer Case – Manufacturing using F<sub>2</sub>

Fluorine substitution in agrochemicals & pharmaceuticals is increasing in popularity however challenges with the safe & selective use of F<sub>2</sub> means using expensive alternative reagents. Sandford & Pichon [1] recently demonstrated synthesis of Flucytosine.

### Challenges in Batch:


- Hazardous reagent (F<sub>2</sub>)
- Highly exothermic reactions (DT up to 180 °C)
- High dilution & long dosing times
- Challenges to perform at scale

### Advantages in Flow:

- Thermal control = Intensification
- Metal-free reactors = Reduced corrosion risk
- Selective mono-fluorination = Increased product purity
- Single reaction step = Reduced isolation costs



A range of products for various application spaces performed in Plantrix®;

- Typically, 10 % F<sub>2</sub> in N<sub>2</sub> & active in polar solvent
- Throughput ranging from 3 to 30 kg/day isolated material (>90 % yield)



## Customer Case – Generic API Manufacturing

Development of a CM route for a marketed drug substance comprising of multiple reactions, crystallisation, filtration & drying!




Scaled-up from Protrix to Plantrix MR555 with ~100 MT/annum capabilities, achieving:

- Process simplification, reduced processing time (52 h to 14 h), improved robustness, reduced cost of goods, power consumption & space utilisation
- Illustrating that generic API's manufactured in Asia still have a lot to benefit from 'going flow'

## Future of Flow Chemistry & Process Intensification

- Supply chain management is favouring in-house / local / distributed manufacturing models
  - Reducing the environmental burden caused by moving raw materials / intermediates / products between sites
- Continuous manufacturing brings an opportunity to further disrupt conventional supply chains, enabling;
  - On-demand manufacturing
  - Reduced 'work in progress' goods
  - Geographic diversification



© Copyright – TNO

- When combined with PAT/model-based predictive control, modular & flexible CM can robustly address the Industry needs!

## Widely Accepted Benefits for Primary Processing

The rationale for a Company to 'Go Flow' is varied & depends on the sector, process type & scale of operation!

### 1. Safe Use of Extreme Reaction Conditions

- Efficient mixing
- Excellent thermal control
- Process intensification of hazardous reactions

### 2. Reduced Development Time

- Small hold-up volume
- Rapid reaction optimisation
- Minimal scale-up steps

### 3. Improved Process Control

- High level of reaction control
- Process reproducibility
- Quality by Design

### 4. Reduced Production Costs

- Increased product quality
- Reduced safety investments
- Higher unit productivity

Engineering for efficiency, flexibility & sustainability bring opportunities to commercialise intensified processes!

- Dramatic reductions in plant size
- Access to reaction conditions previously unobtainable
- Flexibility towards different product types & production volumes




## Where Chemists & Engineers Meet!

“...The only disadvantage of continuous processing, with its different skill & knowledge requirements, is that companies will need to change their ways of working to take full advantage of its benefits.

With typically, 90% of manufacturing costs locked in at the design stage, & the cost to change a process increasing exponentially the further down the development route it is implemented, good management of early process design & development are essential....”



‘The Reality of Continuous Processing’, Manufacturing Chemist, April 2005  
Huw Thomas, Foster Wheeler



Partner to Accelerate Success!

Dr Charlotte Wiles (CEO)  
Chemtrix BV – Headquarters  
Galvaniweg 8A  
6101 XH Echt  
The Netherlands  
E-mail: [c.wiles@chemtrix.com](mailto:c.wiles@chemtrix.com)  
Tel: +31 (0)467 022 600



Please find details of our publications, application notes and white papers on our website:

[www.chemtrix.com](http://www.chemtrix.com)